


dosasap™

A COMPACT TOOL FOR IQ/OQ/PQ PURPOSES

SOFTWARE PACKAGE MAIN FEATURES

- **Dosimetry system calibration** (readout of calibration CTA strips, curve fitting, residual and uncertainty calculations)
- **Electron Beam** width, scan length and scanning uniformity assessment
- Photon and Electron Beam absorbed **dose profile measurement** for dose uniformity estimation in user selected Region Of Interests (ROI) strip length from 5 cm to 12 m

PC controlled dosimetry equipment for CTA dosimeter readout



- **Electron Beam energy determination** with wedge technique according to ISO/ASTM 51649

aerasap™

Software for continuous optical dosimetry tools

data base driven compliant with relevant standards ISO/ASTM 51649, 51650, 51261, ISO 11137-3, ... and FDA CFR 21 part 11.

Contact us

- Sébastien Riegler - Product Manager - Sales - After Sales
s.riegler@aerial-crt.com
 - Abbas Nasreddine - Project Manager - Development
a.nasreddine@aerial-crt.com
- +33 (0)3 88 19 15 15
Parc d'innovation - 250 rue Laurent Fries 67 400 Illkirch - France
www.aerial-crt.com

Aerial
La bonne dose d'innovation
* the best dose of innovation